



Easy to operate with dedicated app

Gel/Polarized light switching

Adopts high-resolution smartphone

Integrated dermatology camera

Handy Dermo

HD-2 : Diagnostic observation required for dermascope possible

- 1) Gel/oil + pressure observation with glass
- 2) Pressure observation with polarizing filter + glass
- 3) Polarizing filter pressureless observation

* Main camera 50 million pixels

* LCD resolution 6.5 inch 1600x720p, Aspect ratio 20:9.

* Standard 10x images can be obtained with a single touch using a dedicated app. Equipped with HDR function to enhance contrast.

* Necessary functions for telemedicine: Wi-Fi video conferencing, telephone use and electric scale.

* Boundary observation and optional normal imaging possible.

◇ Skin contact with silicon adapter
Can be sterilized by autoclaving.

◇ HD-2 Handy Dermo

- 1) Smartphone : Moto G13 (Wide-angle camera: 50 million pixels)
- 2) Size : 186 x 87 x 40.5mm (Including silicon adapter)
- 3) Weight : Single unit: 183.45g, Dermoscope : 224g
- 4) Memory : RAM 4GB+ROM 128GB, micro SD: Max.512GB
- 5) Battery : 5000mAh

◇ Price (All-in-one including smartphone, silicone adapter, light guide, glass, USB cable)

| | |
|------|-------------------|
| | Motorola Moto G13 |
| HD-2 | \$1,000. |

2-11-26 Nagataminami, Minami-ku, Yokohama
Kanagawa, Japan 232-0073

Tel :+81-45-731-2584 Fax: +81-45-714 -3763

E-mail: gawjv440@ybb.ne.jp URL:<http://dermamedical.jp>

Derma Medical Inc.

